

## Cross-Reference-list

Order / article number:

Content/ Volume centrifugetubes	BEKOlut GmbH & Co. KG BEKOlut® QuEChERS	Agilent Technologies Bond Elut®	Sigma Aldrich Supel™ Que	Restek Q-Sep®
<b>EXTRACTION SALTS</b>				
4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g TSCD, 0.5 g DHS	Citrat-Kit-01, <b>CK-01</b> incl. 50mL centrifugetubes	5982-5650		26235
4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g TSCD, 0.5 g DHS	Citrat-Kit-01, <b>CK-01-050</b> without 50mL centrifugetubes	5982-6650	55227-U	26236
4 g MgSO <sub>4</sub> , 1 g NaCl	SALT-KIT-AC, <b>SK-AC</b> incl. 50mL centrifuge tubes	5982-5550		23991
4 g MgSO <sub>4</sub> , 1 g NaCl	SALT-KIT-AC, <b>SK-AC-050</b> without 50mL centrifuge tubes	5982-6550	55294-U	23992
6 g MgSO <sub>4</sub> , 1.5 g NaOAc	SALT-KIT-AC2, <b>SK-AC2</b> incl. 50mL centrifuge tubes	5982-5755		26237
6 g MgSO <sub>4</sub> , 1.5 g NaOAc	SALT-KIT-AC2, <b>SK-AC2-050</b> without 50mL centrifuge tubes	5982-6755	55234-U	26238
<b>DISPERSIVE KITS</b>				
25 mg PSA, 150 mg MgSO <sub>4</sub>	PSA-KIT-01, <b>PK-01</b> in 2mL centrifugetubes	5982-5021		26215
50 mg PSA, 150 mg MgSO <sub>4</sub>	PSA-KIT-01A, <b>PK-01A</b> in 2mL centrifugetubes	5982-5022	55287-U	26124
150 mg PSA, 900 mg MgSO <sub>4</sub>	PSA-KIT-02, <b>PK-02</b> in 15mL centrifugetubes	5982-5056	55228-U	26223
400 mg PSA, 1200 mg MgSO <sub>4</sub>	PSA-KIT-02A, <b>PK-02A</b> in 15mL centrifugetubes	5982-5058	55282-U	26220
25 mg PSA, 25 mg C18e ,150 mg MgSO <sub>4</sub>	PSA-KIT-03, <b>PK-03</b> in 2mL centrifugetubes	5982-5121		26216
50 mg PSA, 50 mg C18e ,150 mg MgSO <sub>4</sub>	PSA-KIT-03A, <b>PK-03A</b> in 2mL centrifugetubes	5982-5122	55288-U	26125
50 mg C18e ,150 mg MgSO <sub>4</sub>	PSA-KIT-03SP, <b>PK-03SP</b> in 2mL centrifugetubes			26242
150 mg PSA, 150 mg C18e ,900 mg MgSO <sub>4</sub>	PSA-KIT-04, <b>PK-04</b> in 15mL centrifugetubes	5982-5156	55229-U	26226
400 mg PSA, 400 mg C18e ,1200 mg MgSO <sub>4</sub>	PSA-KIT-04A, <b>PK-04A</b> in 15mL centrifugetubes	5982-5158	55283-U	26221
25 mg PSA, 2.5 mg GCB ,150 mg MgSO <sub>4</sub>	PSA-KIT-05, <b>PK-05</b> in 2mL centrifugetubes	5982-5221		26217
150 mg PSA, 15 mg GCB ,900 mg MgSO <sub>4</sub>	PSA-KIT-06, <b>PK-06</b> in 15mL centrifugetubes	5982-5226	55230-U	26224
25 mg PSA, 7.5 mg GCB ,150 mg MgSO <sub>4</sub>	PSA-KIT-07, <b>PK-07</b> in 2mL centrifugetubes	5982-5321		26218
150 mg PSA, 45 mg GCB ,900 mg MgSO <sub>4</sub>	PSA-KIT-08, <b>PK-08</b> in 15mL centrifugetubes	5982-5356	55233-U	26225

